Total Hardness Test Kits

The hardness level as mg/L (ppm) calcium carbonate is determined by an EDTA (ethylene-diamine-tetraacetic acid) titration.

HI3812

HI3840

Low Range

HI3841

Medium Range

HI3842

High Range

HI38033

0-30 gpg range



Hydrogen Peroxide Test Kit

The Hanna test kit can quickly and easily determine concentration in water up to 10 ppm of hydrogen peroxide. This method is not affected by stabilizers that can be added to commercial hydrogen peroxide solutions. Hydrogen peroxide concentration is determined by a titrimetric method. It reacts with iodide in acid solution, the amount of iodine generated is equivalent to the hydrogen peroxide in the sample. The liberated iodine is then titrated with standard sodium thiosulfate solution.



Method	Range	Smallest Increment	Chemical Method	# Tests	
HI3812	Hardness (as CaCO ₃ *) Total				
titration	0.0-30.0 mg/L (ppm) 0-300 mg/L (ppm)	0.3 mg/L (ppm) 3 mg/L (ppm)	EDTA	100	
HI3840	Hardness (as CaCO₃*) Total				
titration	0-150 mg/L (ppm)	5 mg/L (ppm)	EDTA	50 avg.	
HI3841	Hardness (as CaCO ₃ *) Total				
titration	40-500 mg/L (ppm)	20 mg/L (ppm)	EDTA	50 avg.	
HI3842	Hardness (as CaCO ₃ *) Total				
titration	400-3000 mg/L (ppm)	100 mg/L (ppm)	EDTA	50 avg.	
HI38033	Hardness (as CaCO₃*) Total				
titration	0-30 gpg	1 gpg	EDTA	100	
HI3844	Hydrogen Peroxide (as H ₂ O ₂)				
titration	0.00-2.00 mg/L (ppm) 0.0-10.0 mg/L (ppm)	0.25 mg/L (ppm) 1.0 mg/L (ppm)	iodometric	100 avg.	
		HI3812 test kit comes with 30 mL hardness buffer, 10 mL calmagite indicator, 120 mL EDTA solution, 20 mL plastic beaker with cap, 50 mL plastic beaker with cap and 1 mL syringe with tip.			
	HI3840 test kit comes with 30 mL hardness LR reagent and 50 mL calibrated vessel.				
	HI3841 test kit comes with 30 mL hardness MR reagent and 50 mL calibrated vessel.				
Ordering Information	HI3842 test kit comes with 30 mL hardness HR reagent and 50 mL calibrated vessel. HI38033 test kit comes with 30 mL buffer solution, 10 mL calmagite indicator, 75 mL EDTA				
	solution (2), 20 mL plastic beaker with cap and 1 mL plastic pipette.				
	HI3844 test kit comes with 100 mL hydrogen peroxide reagent A, 17 g hydrogen peroxide reagent B, 30 mL hydrogen peroxide reagent C, 25 mL hydrogen peroxide reagent D, graduated plastic test tube with cap, 50 mL calibrated plastic vessel, 3 mL plastic pipette, 1 mL plastic pipette and plastic spoon.				
* 1 gpg = 17 ppm CaCO ₂					

