

Specifications	HI701
Range	0.00 to 2.50 ppm
Resolution	0,01 ppm
Accuracy @ 25°C/77°F	±0.03 ppm ±3% of reading
Light Source	LED @ 525 nm
Light Detector	silicon photocell
Environment	0 to 50°C (32 to 122°F); RH max 95% non-condensing
Battery Type	(1) 1.5V AAA
Auto-off	after two minutes of non-use
Dimensions	81.5 x 61 x 37.5 mm (3.2 x 2.4 x 1.5")
Weight	64 g (2.25 oz.)
Method	adaptation of USEPA method 330.5, DPD method
Ordering Information	HI701 Checker®HC is supplied with sample cuvettes with caps (2), free chlorine reagent starter kit (reagents for 6 tests), battery, instructions and quick start guide.

See a list of Checker® reagents and accessories on page 1.24

Free Chlorine

Handheld Colorimeter

- Easier to use and more accurate than chemical test kits
 - EPA approved DPD method
 - ±0.03 ppm ±3% of reading accuracy
 - 0.01 ppm resolution
 - Large, easy-to-read digits
 - Auto shut off

• Dedicated to a single parameter

- Designed to work with Hanna's powder reagents
- Small size, big convenience
 - The Checker®HC easily fits into the palm of your hand or pocket
 - Use for quick and accurate on-thespot analysis
 - Single-button operation: zero and measure

• Ideal for:

- Swimming pools and spas
- Fruit and vegetable sanitation
- Disinfection
- Drinking water and quality control checks

The HI701 Checker®HC bridges the gap between simple chemical test kits and professional instrumentation. Chemical test kits are not very accurate and only give 5 to 10 points resolution, while professional instrumentation can cost hundreds of dollars and can be time consuming to calibrate and maintain. The Hanna HI701 Checker®HC is accurate and affordable.

This Checker®HC portable handheld colorimeter features a resolution of 0.01 ppm and ±0.03 ppm ±3% of reading accuracy. It also uses an EPA approved DPD method.

The contoured style of the Checker®HC fits in your palm and pocket perfectly and the large LCD is easy to read. The auto shut-off feature assures the battery life will not be drained if you forget to turn it off.

The HI701 Checker®HC is extremely simple to use. First, zero the instrument with your water sample. Next, add the reagent. Last, place the vial into the Checker®HC, press the button and read the results. It's that easy.

